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## **AREC**

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# Operating Manual

## **AREC**

### **Heating Magnetic Stirrer with Ceramic Plate**

**F20500010 (230V - 50÷60Hz)**

**F20510010 (115V - 60Hz)**

## **Preface:**

Velp Scientifica thanks you for choosing the heating magnetic stirrer with ceramic plate **AREC**.

The **AREC** model is designed for general use and for all applications which require a set stirring speed and heating plate temperature.

## **Safety instructions:**

Before using the unit, please read the operating manual supplied with the apparatus carefully.

Warning ! Hot surface

Do not dispose of this equipment as urban waste

In order to prevent possible risk of electric shocks, fire and personal injury when the unit is being used, basic safety measures must always be taken, including:

- The heating plate can reach a temperature of up to 540°C. Care must be taken during the heating phase and also when the appliance is switched off as it takes time to cool down.
- The magnetic fields are active on magnetic and metallic parts (magnetic supports, CDs, floppy discs, cardiac stimulators, magnetic cards). Keep these items away from the magnetic stirrer plate and the magnetic stir bars.
- The stirrer must not be used:
  - in hazard areas
  - in bain-marie
  - to stir combustible liquids that have a low combustion temperature
- The vessel must be smaller in diameter than the ceramic plate to prevent any risk of tipping over.
- If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature – the use of Pyrex vessels is recommended.
- The following precautions must be taken before using the instrument:
  - position the instrument on a level surface.
  - position the instrument at least 30 cm from the wall
  - an operator must be present at all times when the instrument is in use
- The liquid must not come into contact with the electric power cable or with the electrical parts inside the instrument.
- Check that the power supply corresponds to the rating plate on the rear of the unit and distant from the heating plate.

- Only correctly earthed three-terminal power cables must be used.
- Do not use the unit if it is not working correctly. In case of malfunctioning, contact your nearest assistance centre.
- Personal protective equipment must be compatible with the possible risks posed by the material being processed and the glass parts.
- Follow the cleaning instructions described in this manual

**This unit must only be used for laboratory applications**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

This unit has been designed and produced in compliance with the following standards.

Safety requirements for electrical apparatus for:

Measurement and control and for laboratory use	CEI EN 61010-1
Electrical equipment for laboratory use	UL 3101-1
General requirement – Canadian electrical code	CAN/CSA-C22.2

***N.B.:*** *The manufacturer is committed to constantly improving the quality of the products and reserves the right to modify the characteristics without prior notice.*

# CONTENTS

<b>1</b>	<b>INTRODUCTION</b>	<b>1</b>
1.1	Parts included	2
1.2	Instrument description	2
1.3	Knobs and functions	3
1.4	Information regarding construction materials	3
<b>2</b>	<b>ASSEMBLY AND INSTALLATION</b>	<b>4</b>
2.1	Electrical connection	4
2.2	Start-up	4
<b>3</b>	<b>OPERATING CONTROLS</b>	<b>4</b>
<b>4</b>	<b>END OF WORK CYCLE OPERATIONS</b>	<b>4</b>
<b>5</b>	<b>MAINTENANCE</b>	<b>5</b>
5.1	Cleaning	5
5.2	Disposal	5
5.3	Error messages	6
<b>6</b>	<b>ACCESSORIES AVAILABLE ON REQUEST</b>	<b>7</b>
<b>7</b>	<b>SPARE PARTS</b>	<b>8</b>
<b>8</b>	<b>TECHNICAL FEATURES</b>	<b>8</b>
<b>9</b>	<b>WIRING DIAGRAM</b>	<b>9</b>
<b>10</b>	<b>WARRANTY</b>	<b>10</b>

---

# 1. Introduction

---

The **AREC** heating magnetic stirrer with ceramic plate, is used in the laboratory for general use and for all the applications that require a precise regulation of the stirrer speed and the heating plate temperature.


The ceramic heating plate **(1)** makes it possible to obtain higher temperatures than with the traditional aluminium heating plate and a faster heating of the sample. Ceramic is an inert and very hard material and is resistant to almost any type of chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The AREC is used to mix liquids using the magnet drive of a stir bar placed inside the vessel. The size and shape of the magnetic stir bar determines the efficiency of the stirring itself at any given stirring speed. It is therefore advisable to choose the most suitable magnetic stir bar in relation to the quantity and characteristics of the liquid and the type of vessel used.

The stir bar which satisfies most applications is code A00001060 (10x40 mm).

The stirring speed **(2)** and the temperature **(3)** of heating plate, indicated on the digital display **(4)**, are adjusted using the relative knobs on the frontal panel. The “Hot Plate” light **(5)** turns on when the temperature of the heating plate is higher than 50°C.

**NOTE:** If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

## 1.1 Parts included

Check the unit is complete after unpacking.  
The table below shows the parts included.

Parts included		Qty
1	Heating magnetic stirrer	1
2	Power cable	1
3	Operating manual	1

## 1.2 Instrument description

The structure of the unit, made out of die-cast aluminium, combined with its innovative design and low profile, allows the operator to work in extreme comfort and safety under all conditions.

The new design combines ergonomics with a high Electrical Protection rating (IP42) according to the rules CEI-EN 60529.

The control panel that contains the electronic parts, is completely separate from the heating plate in order to increase the safety rating during use. The inclination of the knobs has been carefully studied in order to facilitate use by operator.

The stirring speed can be set within a range of 50-1300 rpm using the knob on the front panel.

The AREC has an electronic feed-back speed control to ensure a constant stirring speed even if there are variations in the working conditions (i.e. sample viscosity, friction of stirring bar inside the container, etc.).

The temperature of the heating plate can be set within a range of 5 – 540°C and with a 5°C resolution, using the knob on the front panel.

The temperature setting is shown on the display.

A temperature probe in direct contact with the heating plate takes continuous temperature readings whilst an electronic closed-loop feedback system ensures that the set temperature is maintained.

## 1.3 Knobs and Functions

The knobs are positioned on the front panel, are easily accessible and well away from the heating plate in order to ensure maximum operator safety as well as safeguarding the electronic components inside the unit.

<p><b>STIRRING SPEED</b></p> 	<p>Turn the stirring speed knob/switch on the front panel in order to start the stirring function. It is possible to select speeds from between 50 and 1.300 rpm.</p> <p>If the heating function is off, the display indicates that stirring is operational by bright bars rotating on the display. The AREC has an electronic feed-back speed control to ensure a constant stirring speed even if there are variations in the working conditions.</p>
<p><b>TEMPERATURE</b></p> 	<p>Turn the heating knob/switch on the front panel in order to start the heating function. It is possible to select heating plate temperatures from 5 to 540°C with a 5°C resolution. The set temperature is indicated on the display.</p>
<p><b>“HOT PLATE”LIGHT</b></p> 	<p>The HOT PLATE light is always on when the temperature of the heating plate exceeds 50°C.</p> <p>When the heating function is turned Off (Heating knob on OFF) and the heating plate temperature is over 50°C, the HOT PLATE light flashes and OFF is indicated on the display until the temperature of the heating plate falls below 50°C.</p> <p><b><i>N.B.</i></b> This warning is not active if the instrument is not connected to the electric power supply..</p>

## 1.4 Information regarding construction materials

Housing	Die-cast aluminium
Paint	PET
Front control panel	PET
Feet	EPDM
Heating plate	Ceramic

If any liquid being processed falls onto the above materials, clean them immediately to avoid chemical corrosion.

---

## 2. Assembly and installation

---

### 2.1 Electrical connection

After having unpacked the instrument, place the unit on a laboratory bench. Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the cable supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

### 2.2 Start-up

After the connecting the instrument to the power supply, place a flask containing the sample and a suitable magnetic stir bar on the heating plate.

---

## 3. Operating controls

---

Start the stirring function by turning the stirrer knob “Stirrer rpm” on the front panel. Speeds from 50 to 1.300 rpm can be selected using the analogical scale around the knob.

Bright bars rotating on the display indicate that the instrument is running. The speed of rotation of the bright bars on the display is proportional to the stirring speed.

*NOTE: The rotating bars are not display when the temperature function is in use.*

To start the heating function turn the heating knob “Heating °C” on the front panel. Temperatures from 5 to 540°C can be selected at 5° intervals. The temperature setting is shown on the display.

The HOT PLATE light remains lit when the temperature of the heating plate exceeds 50°C.

When heating is turned Off (Heating knob OFF) and the heating plate temperature is over 50°C, the HOT PLATE indicator flashes and OFF is shown on the display until the temperature of the heating plate falls to below 50°C.

**N.B.** *This warning is not active if the instrument is not connected to the power supply.*

---

## **4. End of work cycle operations**

---

At the end of the work cycle turn off the stirring and heating knobs.

---

## **5. Maintenance**

---

No routine or extraordinary maintenance is necessary apart from periodically cleaning of the unit as described in this manual (chap.5.1).

In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

Should the fuse need changing, disconnect the unit from the power supply and open the cover on the back of the unit to access the two fuses 5A (230V) or 8A (115V).

### **5.1 Cleaning**

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable and non-aggressive detergent to clean the unit.

### **5.2 Disposal**

The final disposal of the unit or of its components must be carried out in compliance with the directives and laws in force in the country concerned.

## 5.3 Error messages

If the display shows an error message please contact Velp Scientifica's technical service department.

When the display shows an error message, the stirring and heating functions stop automatically.

**NOTE:** To remove the error message disconnect the instrument from the power supply.

Message	Problem
AL1	Thermocouple is out of range
AL2	Excessive heating time
AL3	The stirring system doesn't run correctly

## 6. Accessories available on request

Description:	Code No
Magnetic stir bar Ø. 6x35mm	A00001056
Magnetic stir bar Ø 10x40mm*	A00001060
Magnetic stir bar Ø 9.5x60mm	A00001061
Hemispheric bowl for 250ml flasks	A00001071
Hemispheric bowl for 500ml flasks	A00001072
Hemispheric bowl for 1000ml flasks	A00001073
Support rod	A00001069

\* Magnetic stir bar recommended for AREC

		
A00001056	A00001060	A00001061
		
A00001071	A00001072	A00001073
		
A00001069		


## 7. Spare parts

Description:	Code No.
Knob	10002097
Foot	10000239

## 8. Technical features

Power supply	V/Hz	230/50-60 115/60
Total electric power	W	800
Dimensions	mm (wxhxd)	205x96x335
	Inch (wxhxd)	8x3.7x13.2
Weight	Kg	3.9
	lbs	9.6
Speed range	rpm	50-1,300
Speed control		Electronic feedback
Speed selection		Analogue
Temperature range	°C	5-540
	°F	41-1004
Temperature control		Electronic probe
Temperature selection		Display reading 5°C intervals
Max load-bearing capacity of plate	Kg	25
	lbs	67
Ceramic plate dimensions	mm	200x200
Construction material		Die cast aluminium
Environmental temp. range	°C	+5...+40
	°F	+41...+104
Storage temp. range	°C	-10...+60
	°F	+14...+140
Max. humidity	%	Max 80
Operating mode		Continuous
Electrical safety level CEI EN60529	I	P42
Pollution degree CEI EN61010-1		2
Installation category CEI EN61010-1		2
Max. altitude	m	2000
Fuse	A	2x5A (2x8A – 115V/60Hz)

## 9. Wiring diagram

- 1) Main board
- 2) Thermocouple
- 3) Stirrer potentiometer
- 4) Temperature potentiometer
- 5) Display board
- 6) Revolution counter board
- 7) Heating element
- 8) Motor
- 9) Fuses rit. 5x20mm n.2x5A (n.2x8A-115V)

---

## 10 Warranty

---

The unit is guaranteed against production defects for 25 months from our invoice date.

In accordance with this guarantee VELP SCIENTIFICA undertakes to repair any units resulting as faulty due to the quality of the materials used or poor workmanship.

Units rendered faulty due to inexpert handling/use or carelessness will not be replaced or repaired under warranty.

For more details please contact your Distributor.

### **Exclusions:**

The guarantee will be considered null and void for faults resulting from:

- inexperience and carelessness of the operator
- repairs, maintenance or replacement of parts carried out by personnel or Companies not authorized by the manufacturer
- use of the instrument that does not comply to the instructions/recommendations given in the present operating manual
- use of non-original spare parts.

## Declaration of conformity

We the manufacturer	VELP SCIENTIFICA s.r.l.
Address	Via Stazione, 16 USMATE (Milan) Italy

under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1 (2001)

EN 61326-1 (1997) + A1 (1998) + A2 (2001) + A3 (2003)

2002/95/EC (RoHS)

2002/96/EC (WEEE)

and satisfies the essential requirements of the following directives:

Machinery directive 2006/42/EC

Low voltage directive 2006/95/EC

Electromagnetic compatibility directive 2004/108/EC

plus modifications and that the documents listed in annex I are available at Velp's offices as foreseen by the machinery directive.

## ***Thank you for having chosen a VELP product !***

*Since 1983 Velp has offered to professionals in the sector a range of sophisticated and reliable equipment with excellent operating capacity and backed by high levels of know-how at competitive prices.*

*Velp works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification. Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.*

*Let us present our product Lines:*

### ***Environment***

*THERMOREACTORS ECO SERIES  
B.O.D. DETERMINATION SYSTEMS  
REFRIGERATED THERMOSTATS  
COOLED INCUBATORS  
FLOCCULATORS  
OVERHEAD MIXER  
MINERALIZATION SYSTEM FOR HEAVY METALS TRACE DETERMINATION  
TURBIDIMETER  
RADIATION DETECTOR*

### ***Stirring***

*HEATING MAGNETIC STIRRERS  
VERTEX DIGITAL THERMOREGULATOR  
ULTRAFLAT MAGNETIC STIRRER  
MAGNETIC STIRRERS  
OVERHEAD STIRRERS  
HEATING PLATES  
VORTEX MIXERS  
HOMOGENIZER*

### ***Food&Feed***

*DIGESTERS DK SERIES  
JP RECIRCULATING WATER PUMP FOR FUMES ASPIRATION  
SMS SCRUBBER  
DISTILLATION UNIT UDK SERIES  
SOLVENT EXTRACTORS  
RAW FIBER EXTRACTORS  
DIETARY FIBER EXTRACTORS  
LABORATORY AIDS  
OXITEST OXIDATION TEST REACTOR*

### ***Pumps***

*RECIRCULATING WATER VACUUM PUMP  
PERISTALTIC PUMPS*

**Fisher Scientific**

Part of Thermo Fisher Scientific

<b>Belgique / België</b>	tél 056 260 260	fax 056 260 270	be.fisher@thermofisher.com	www.be.fishersci.com
<b>España</b>	tfno 91 380 67 10	fax 91 380 85 02	es.fisher@thermofisher.com	www.es.fishersci.com
<b>France</b>	tél 03 88 67 53 20	fax 03 88 67 11 68	fr.commande@thermofisher.com	www.fr.fishersci.com
<b>Italia</b>	tel 02 953 28 258	fax 02 953 27 374	it.fisher@thermofisher.com	www.it.fishersci.com
<b>Netherlands</b>	tel 020 487 70 00	fax 020 487 70 70	nl.info@thermofisher.com	www.nl.fishersci.com
<b>Portugal</b>	tel 21 425 33 50	fax 21 425 33 51	pt.fisher@thermofisher.com	www.pt.fishersci.com
<b>Suisse / Switzerland</b>	tél 056 618 41 11	fax 056 618 41 41	info.ch@thermofisher.com	www.ch.fishersci.com

**Velp Scientifica srl**  
20040 Usmate (Milan) Italy  
Via Stazione, 16  
tel. +39 039 628811  
fax +39 039 6288120  
e-mail: [inse@velp.it](mailto:inse@velp.it)  
Internet: [www.velp.com](http://www.velp.com)

Distributed by: